

FA Stand for Zeiss Stemi 2000

The FA Stand for Zeiss Stemi 2000 was developed for the inspection of larger samples, which may not be placed on standard microscope stands. For microscopic analysis the samples lie directly on the table. The stand offers maximum stability and flexibility and is ideally suited for use in routine and research laboratories.

**FA Stand with
Zeiss Stemi 2000 C**



Description:

- Base plate, stain, painted black, dimensions: 250 x 250 x 25 mm
- Vertical column long, \varnothing 32 mm, length 500 mm
- Spacer 100 mm to extend the distance between optical axis and column, with star grip for clamping
- Vertical column short, \varnothing 32 mm, length 100 mm, to adapt Stemi mount
- Security cap for vertical column

Order Information:

Art.-No.. 3055281

FA Stand FA5281 for Zeiss Stemi 2000

OptoSys GmbH

Technologie- und Innovationszentrum
Robert-Bosch-Straße 7
64293 Darmstadt

Tel. +49 (0)6151 - 44 329
Fax +49 (0)6151 - 42 66 43
Email: info@opto-sys.de
www.opto-sys.de

Optische Komponenten &
Systemlösungen GmbH

