



Universal Stand with Sliding Stage for Zeiss SteREO Discovery and Axio Zoom.V16

The Universal Stands for Zeiss SteREO Discovery and Axio Zomm.V16 have been developed for the inspection of larger samples, which may not be placed on standard microscope stands. This Universal Stand for reflected illumination is equipped with a manual x/y sliding table with a large travel range.



Zeiss SteREO Discovery.V20 and Universal Stand with Sliding Table

Description:

- Base plate and sliding table aluminium, black anodized
- Dimensions of base plate: about 600 x 400 x 25 mm
- Travel range about 300 x 300 mm
- Extension of the distance between focusing device and optical axis by installing a 150 mm spacer
- A motorized focusing drive is recommended

Order Information:

Art. No. 3055262 Universal Stand A5262 with sliding table

OptoSys GmbH

Technologie- und Innovationszentrum
Robert-Bosch-Straße 7
64293 Darmstadt

Tel. +49 (0)6151 - 44 329
Fax +49 (0)6151 - 42 66 43
Email: info@optosys.eu
www.optosys.eu www.kursmikroskope.eu

Optische Komponenten &
Systemlösungen GmbH

