

EasyLED Spotlight Plus

Ergonomic and compact illumination system for stereo microscopy and macroscopy



The SCHOTT **EasyLED Spotlight Plus** System has been specifically developed for microscopy illumination tasks. Each lighthouse offers powerful 140 lm and due to its special construction allows optimal heat dissipation over the cooling fins so that the illuminated objects will not be heated. For the Double Spotlight Plus the user can select one or both lighthouses for ideal contrast by repeatedly pressing the On-button. For the integrated ergonomic control unit various mounting adaptors and stand-alone solutions are available.

The EasyLED Spotlight Plus Systems comes with an external wide range of power supply (100-240 V; 50-60 Hz). It has been approved and certified as laboratory equipment and is in compliance with the CE regulations.

Features:

- Spotlight Plus adaptors for all stands plus stand-alone version available
- Double Spotlight Plus with switching feature: LED head 1 ON, LED head 2 A, both LED heads ON, both LED heads off, etc.
- Max. light flux: 130 lm each spot